


RENAL
CARE



D.med NephroFlow™

RECIRCULATION AND ACCESS
FLOW MEASURING DEVICE



BETTER LIFE EXPECTANCY TO THE VASCULAR ACCESS OF YOUR PATIENT



LARGE AREA TOUCH SCREEN

The high-sensitivity screen optimizes usability during the treatment. Intuitive interface guides the user through measurement and provides all the steps and advices.



QUICK MEASUREMENT MODE

The quick mode allows the user to perform a measurement without entering patient data. D.med NephroFlow is ready to measure within seconds.



BUILT-IN FAST START

It takes only few seconds to be ready to measure once the device is switched on. No loss of time. Immediate effectiveness.





6 MONTHS TREND VIEWER

Graphical documentation of previous measurements up to 6 months allows to supervise the vascular access of the patient and helps for proactive minimal intervention.



PRINT & SAVE VIA USB

In Normal mode, all the data is immediately stored in the patient record. A simple USB connection allows you to export and/ or print the data without use of an external PC.



Vascular dialysis access is considered to be the Achilles heel of hemodialysis. Following the advises of the NKF KDOQI™, a surveillance to diagnose and prevent the progression of stenosis or thrombosis is recommended at least once a month.*

D.med NephroFlow permits the needed surveillance, using the ultrasound dilution technology, a noninvasive access flow monitoring. Safe, quick and accurate, the device provides the actual blood flow, the access recirculation and the shunt flow of your patient during hemodialysis.

Due to its high quality engineering, on D.med NephroFlow the annual functional inspection can easily be done on site, without the need to remove the device from its operating environment. Furthermore, our tubing database is constantly updated and can be simply downloaded when needed.

*National Kidney Foundation: KDOQI clinical practice guidelines for vascular access. Updated 2006



Nipro Medical Europe NV
European Headquarters, Weihoek 3H, 1930 Zaventem, BELGIUM
T: +32 (0)2 725 55 33 | F: +32 (0)2 725 70 41
info@nipro-mediceurope.com | www.nipro-mediceurope.com

NIPRO EUROPE - FRANCE : Biopôle Clermont-Limagne, 63360 Saint Beauzire, FRANCE | T: +33 (0)473 33 41 00 | F: +33 (0)473 33 41 09
NIPRO D.MED GERMANY GMBH : Bernhard-Nocht-Str.99, Hamburg 20359, GERMANY | T: +49 (0)40 75 66 32 50 | F: +49 (0)40 75 66 325 100
D.MED HEALTHCARE AG LTD. : HNO B1, Sowutuom, PO Box CT 1919 Accra, GHANA
NIPRO EUROPE - ITALY : Centro Direzionale Milanofiori, Strada 1 - Palazzo F/1, 20090 Assago (Milano), ITALY | T: +39 (0)2 57 50 00 57 | F: +39 (0)2 57 51 81 11
NIPRO D.MED AG KAZAKHSTAN : Gogol Street 39, of. 325, 050044 Almaty, KAZAKHSTAN
NIPRO EUROPE - NETHERLANDS : Regus Business Center, Verlegde Poolseweg 16, 4818 CL Breda, NETHERLANDS | T: +31 (0)76 524 50 99 | F: +31 (0)76 524 46 66
NIPRO D.MED NIGERIA LTD. : 5 Omowunmi Street, Oke-Afa, Isolo, Lagos, NIGERIA
NIPRO EUROPE - POLAND : Ul. Panska 73, 00-834 Warszawa, POLAND | T: +48 (0)22 31 47 155 | F: +48 (0)22 31 47 152
NIPRO EUROPE - RUSSIA : 12 Krasnopresnenskaya Nab., Office 1408, entrance 6, 123610 Moscow, RUSSIA | T: +7 (0)495 258 1364 | F: +7 (0)495 258 1367
NIPRO SENEGAL S.U.A.R.L. : 27 Avenue Georges Pompidou, Dakar, SENEGAL
NIPRO D.MED D.O.O. BELGRADE : Bastovanska 68a, 11000 Belgrade, SERBIA | T: +381 (0)11 75 15 578
NIPRO D.MED HEALTHCARE SOUTH AFRICA (PTY) LTD : Unit 2, 4B Dwyka Street, Stikland Industria, Cape Town, 7530, SOUTH AFRICA | T: +27 21 949-2635 | F: +27 21 949-2397
NIPRO D.MED HEALTHCARE SOUTH AFRICA (PTY) LTD : 226 Market Str, Fairland, Randburg, Johannesburg, 2170, SOUTH AFRICA | T: +27 11 431 1114 / 26 | F: +27 11 431 1115
NIPRO EUROPE - SPAIN : Poligono Los Frailes n° 93 y 94, Daganzo, 28814 Madrid, SPAIN | T: +34 (0)91 884 5531 | F: +34 (0)91 878 2840
NIPRO NORDICS A.B. : Stortorget 13B, 211 22 Malmö, SWEDEN | T: +46 (0)40 60 61 250 | F: +46 (0)40 60 61 255
NIPRO D.MED SCHWEIZ GMBH : Kantonsstrasse 77, 8807 Freienbach, SWITZERLAND | T: +41 (0)55 410 40 44 | F: +41 (0)55 410 40 43
NIPRO D.MED SAGLIK HIZMETLERI TIC. LTD. STI. : Aziziye Mah.Pilot Sok.17/6, Çankaya-ANKARA, TURKEY | T: +90 (0)312 442 21 12 | F: +90 (0)312 442 21 92
NIPRO D.MED UK LTD. : 25 Barnes Wallis Road, Segensworth East, Fareham Hampshire PO15 5TT, UNITED KINGDOM | T: +44 148 985 48 30

Products are continuously under review in the light of technical advancement; the actual specifications may therefore be subject to improvement or modification without notice. Some features and/or accessories displayed and described could be optional on different markets. The information provided does not constitute an offer to sell and is provided for merely informational purposes. Images published in the present document are shown exclusively for demo purposes and are representative of the possible applications and might be adapted for graphic needs. All rights reserved.

Nephroflow-EN-04_May.16