

Hemodialysis is easier now with customized kits

We offer single use sets for connection & disconnection procedures in each hemodialysis fistula and catheter therapy. These sets are made of supplementary components that create complete kits tailored to your individual needs and practices while providing you the safest and the easiest solution.

MAKING YOUR LIFE EASIER WITH A WIDE RANGE OF ON/OFF KITS

Customized series of procedure sets are exceptional time and costs savers not only for medical professionals but also for administrative and support services.

Benefits of using customized kits:

- Reduced preparation times as all items are ready to use at the exact needed time: preparation time decreases approximately 50% compared to traditional preparation methods.
- Reduced risk of cross infection as all needed components are ready to use & sterile in a single package
- The sets include all necessary elements to ensure optimally prepared procedures and avoid stress due to lack of the required items during interventions
- Efficient warehouse management that simplifies the ordering process: one set means one procedure
- Financial and administrative optimization

KITS EXAMPLES

CVC CONNECTION/ DISCONNECTION KIT

Connection content:

- · 2x Nitrile Glove, M
- 1x "T" Split Drape, 30x40 cm ±2 cm, T split part 11 cm ±1 cm
- 4x Alcoholic 2% Chlorhexidine Swab, 6x6 cm nonwoven fabric 9x12 ±1 cm
- 2x Luer Slip Syringe, 5 ml

Disconnection content:

- 2x Nitrile Glove, M
- 4x Alcoholic 2% Chlorhexidine Swab, 6x6 cm nonwoven fabric 9x12 ±1 cm
- 1x Catheter Pouch 18 ±1 cm x 7 cm ±0,3 cm with plaster included, 10x7,5 ±0,3 cm



BASIC CENTRAL VENOUS CATHETER-CVC- INSERTION PACK

Content:

- 1x Self Adhesive Cover Drape, 90x75 cm
- 1x CVC Drape with Adhesive Fen, 100x120 cm
- 1x Scalpel, No 11
- 5x Gauze Swabs, 10 cm x 10 cm
- · 2x Gallipots, 120 ml
- 2x Sponge Stick
- 1x Syringe, 10 ml, Luer Slip
- 1x Hypodermic Needle 21G 38 mm Green
- 1x Mersilk Suture, W775H
- 1x Tray, 32x15x5 cm



FISTULA CONNECTION/ DISCONNECTION KIT

Connection content:

- · 2x Nitrile Glove, M
- 1x Plain Drape, 50x40 ±2 cm
- 2x Alcoholic %2 Chlorhexidine Swab, 6x6 cm, nonwoven fabric 9x12 ±1 cm
- 2x Non-woven Swab,
 7,5x7,5 cm ±0,5 cm, 4 ply
- 4x Plaster tape, 1,25x13 ±0,3 cm Notch part 1 cm ±0,3 cm
- 4x Plaster tape, 2,5x13 ±0,3 cm, Notch part 1 cm ±0.3 cm



Disconnection content:

- · 2x Nitrile Glove, M
- 1x Vinyl Glove, L
- 5x Non-woven Swab,
 7,5x7,5 cm ±0,5 cm, 4 ply
- 2x Wound Pad 5x7,5 cm

DIALYSIS FISTULA CONNECTION/ DISCONNECTION KIT

Connection content:

- 1x Plain Drape, 50x40 cm
- 2x Alcohol Swab, 6x6 cm
- 8x Plaster Tape, 2,5x15 cm

Disconnection content:

- 4x Swab Ball
- 2x Fistula Pressure Bandage



BUILD YOUR OWN PROCEDURE

With a comprehensive set of products that meets all EU standard quality requirements, we enable you to customize and create your own set.

Some of the items you can include in your personalized kits are:

- Nitrile gloves
- Plaster strips
- Field drapes
- Catheter drapes,
 T-split
- Non-woven swabs
- Alcoholic chlorhexidine swab

- Fistula pressure bandages
- Single use tourniquets
- Wound pads
- Masks
- Syringes
- ... And much more.



High quality of tailor-made solutions

Our team makes an extraordinary commitment to deliver **high-quality** and safe solutions in the easiest way possible. We ensure that the kits meet and exceed customers' and patients' and expectations of quality and safety.

The manufacturing process takes place under EU regulations and ISO standard conditions.



Contact us for more information at:



Distributor:

D. Med Healthcare Germany AG Klaus Bungert Str. 3 40468 Dusseldorf, Germany Phone: +49 211 65 04 15 60 Fax: +49 211 65 04 15 61

sales@dmed-healthcare.com www.dmed-healthcare.com